



# Antimicrobial Medical Adhesive tape

A critical technological advance in the field of medical adhesive tape: TRIOMED™ incorporates a powerful antimicrobial to control microbiological contamination (fungi, bacteria, viruses)



## Indications for use:

- The TRIOMED™ Medical Adhesive Tape is a single use, disposable device intended for general securement of tubing, wound dressings, electrodes or any other device attachment.
- Latex free
- Made with a conformable and flexible non-woven
- Provides a strong and reliable fixation, comes in rolls and is easy to cut, remove paper backing and place over intended area
- 5 year shelf life (proven efficacy)
- NO more wastage

## Laboratory tested advantages of the TRIOMED™ Antimicrobial Medical Adhesive Tape

- Proven to kill on the external surface of the tape at least 99.99% of Staphylococcus aureus MRSA, Enterococcus faecalis VRE, Klebsiella Pneumoniae, Pseudomonas aeruginosa, Escherichia coli, Acinetobacter baumannii, and Influenza A H1N1
- Non-cytotoxic and non-irritating
- Releases no chemicals on the patient
- Hypoallergenic



### SIZES:

- 2.5cm x 10m
- 5cm x 10m
- 10cm x 10m



Health Canada / Santé Canada



BIOMEDICAL

141 63 boul. Du Curé Labelle Suite 50  
Mirabel, Québec, Canada J7J 1M3  
| p: 438.792.6288

WWW.I3BIOMEDICAL.COM

## Neutralizes 99.99% of harmful pathogens on contact

All Scientific and clinical studies indicate that the external surfaces of medical & surgical tapes and dressings commonly used in healthcare settings are contaminated with pathogenic bacteria and may serve as a significant source of infection.

The patented TRIOMED™ technology incorporates a broad-spectrum and powerful Tri-iodide antimicrobial engineered to eliminate this infection risk.


The TRIOMED™ Antimicrobial Medical Adhesive tape will effectively kill on its external surface harmful infections and is the only existing solution to stop this widespread contamination.

**Study comparing the microbiological contamination of commercially available medical tapes versus TRIOMED™ tape after only 1 hour on a patient.**



\* Microbiological contamination of Commercially available tapes

VS



\* NO microbiological contamination on TRIOMED™ Antimicrobial Medical Adhesive tape

**Conclusion:** As evidenced by pictures above, the commercially available tape's external surface is microbiologically contaminated and being digested by the infection source, while the TRIOMED™ tape maintains its microbiological integrity.

**REFERENCES:** [Scientific publications on medical tape and dressing contamination:](#)

1. Redelmeier, DA and Livesley, NJ, Adhesive Tape and Intravascular-Catheter Associated Infections. J Gen Intern Med. Vol. 14, p. 373-375, 1999.
2. Lavelle BE. Reducing the Risk of Skin Trauma Related to Medical Adhesives. Managing Infection Control. June 2004.
3. Harris PNA, et al. Adhesive Tape in the Health Care Setting: Another High-Risk Fomite? Medical Journal of Australia. Vol. 196:1, p. 34, Jan. 16, 2012.
4. Berkowitz DM, et al. Adhesive Tape: Potential Sources of Nosocomial Bacteria. Applied Microbiology. Vol. 28, No. 4. P. 651-654, October 1974.
5. Wilcox MH, et al. A Five Year Outbreak of Methicillin-Susceptible Staphylococcus aureus Phage Type 53,85 in a Regional Neonatal Unit. Epidemiol Infect. Vol. 124. P. 37-45, 2000.
6. Dickinson M, et al. Diagnosis and Successful Treatment Complicating Endotracheal Intubation: Cutaneous Zygomycosis (Mucormycosis). Chest. Vol. 114. p. 340-342, 1998.
7. Everett ED, et al. Rhizopus Surgical Wound Infection Associated With Elasticized Adhesive Tape Dressings. Arch Surg. Vol. 114. P. 738-739, 1979.
8. Arias KM. Contamination and Cross Contamination on Hospital Surfaces and Medical Equipment. Initiatives in Safe Patient Care. Accessed at: [www.initiatives-patientsafety.org](http://www.initiatives-patientsafety.org)
9. Cady, M, DO , Gross ,L, Lee, I, V Tape: A potential vector for infection. J.APSF, 2011
10. G.Christiaens ,M.P Hayette, D.Jacquemin ,P.Melin ,J.Mustsers ,P. De Mol : An outbreak of Absidia Corymbifera infection associated with bandages contamination in a burns unit , The Journal of Hospital infection , September 2005, volume 61, issue 1 ,P.88